

TK019

臨床データ

10in1	Flu A	Flu B	RSV	ADV	MP
Sensitivity 感度	98.10% (103/105)	98.91% (91/92)	96.97% (96/99)	98.85% (86/87)	95.61% (109/114)
Specificity 特異度	95.65% (154/161)	95.40% (166/174)	98.11% (156/159)	97.08% (166/171)	96.40% (107/111)
Total Coincidence Rate 全体一致率	96.62% (257/266)	96.62% (257/266)	97.67% (252/258)	97.67% (252/258)	96.00% (216/225)
10in1	hRV	COVID-19	BoV	HPIV	BP
Sensitivity 感度	96.97% (96/99)	95.61% (109/114)	96.63% (86/89)	98.39% (61/62)	96.94% (95/98)
Specificity 特異度	98.08% (153/156)	96.45% (136/141)	96.99% (161/166)	97.41% (188/193)	96.82% (152/157)
Total Coincidence Rate 全体一致率	97.65% (249/255)	96.08% (245/255)	96.86% (247/255)	97.65% (249/255)	96.86% (247/255)