

臨床データ

Candida albicans test results:

	Cultivation method test results		Total
	Positive	Negative	
Positive	118	2	120
Negative	1	229	230
	119	231	350

臨床感度は99.16% (95%CI=95.39%～99.85%)、臨床特異度は99.13% (95%CI=96.90%～99.76%)、臨床合計一致率は99.14% (95%CI=97.51%～99.71%) です。

Trichomonas vaginalis test results :

	Cultivation method test results		Total
	Positive	Negative	
Positive	109	1	110
Negative	2	238	240
	111	239	350

臨床感度は98.20% (95%CI=93.67%～99.50%)、臨床特異度は99.58% (95%CI=97.67%～99.93%)、臨床合計一致率は99.14% (95%CI=97.51%～99.71%) です。

Gardnerella vaginalis test results:

	Cultivation method test results		Total
	Positive	Negative	
Positive	145	1	146
Negative	1	203	204
	146	204	350

臨床感度は99.32% (95%CI=96.22%～99.88%)、臨床特異度は99.51% (95%CI=97.82%～99.91%)、臨床合計一致率は99.43% (95%CI=97.94%～99.84%) です。